

Are your PCR kits
**CUSTOMER
FRIENDLY?**

MATCH
biosystems

MATCH uxPCR

ROOM temperature stability



Dx industry

Currently		INDUSTRY BREAKTHROUGH ROOM temperature Ready-To-Use
Shipping-Storage -20°C		MATCH uxPCR
INEFFICIENT 💰 ⏳ 🤒		
• Dry ice shipping. • Exp. in customs. • No RTU ➔ Technical errors ➔ Contamination ➔ False positives		
PERFORMANCE (sensitivity, LOD)	Decreased	Max
MANUFACTURING times	2-3 days	< 1 hour
Initial INVESTMENT	>100K€ (small batches)	0 € (OEM service)
IN-HOUSE tuning	Strongly hard (expensive)	N/A
Energy consumption	High	Low (SDG)
Chemicals with greenhouse effect	Cooling fluids	FREE (SDG)
		>90% OP COST REDUCTIONS 100 plates/day-shift



matchbios@matchbios.com
+34 626 660 555

MATCH uxPCR

Are your PCR kits customer friendly?

SUSTAINABLE
DEVELOPMENT

GOALS

3

GOOD HEALTH
AND WELL-BEING

9

INDUSTRY, INNOVATION
AND INFRASTRUCTURE

13

CLIMATE
ACTION




Potential in Dx industries

INDUSTRY BREAKTHROUGH

ROOM temperature
Ready-To-Use

MATCH uxPCR




matchbios@matchbios.com
+34 626 660 555

ASK ME FOR A PARTNERSHIP



Adrián H Teruel, PhD
CEO

adrian@matchbios.com
<https://www.linkedin.com/in/adrianhteruel/>
<https://matchbios.com>
+34 626 660 555



CERTIFICATE NUMBER

MD 798249



NÚMERO DE LICENCIA

8937-PS

CSV: ZWC2YPB90C



INNOVATIVE SME



MATCH

biosystems

SUSTAINABILITY THAT SAVES THE PLANET
TIME THAT SAVES LIVES